



# 5-second Effort ... a Perfectly-sealed Petri Dish




## Shrink Seals



### Fast & Easy to Use

- Slip moist band over Petri dish edge
- Each plate is 5 seconds of effort for 100% security
- No wrapping of dish required



### Safe & Secure

- Tight, smooth seal every time
- No gaps or ripping
- Pathogenic microorganisms safely contained



### Preserves Media Integrity

- Blocks evaporation during extended incubation
- Extends shelf-life of sensitive or pre-made plates

*Other products available at  
Scientific Device Laboratory*


#### TIME-SAVING MYCOLOGY PRODUCTS

- *Avoid bubbling or foggy tape preps with Fungi-Tape (#745)*
- *Create beautiful permanent mounts with built-in stain MycoPerm (#753)*
- *Shrink Seals (#261)*

*For more  
information  
see the  
catalog on  
our website.*



# Time-saving Mycology Products that break the mold.



## MycoPerm

- Integrate your permanent mycology mounting & stain into one step
- Protect preps with plasticizers for a durable, long-life with this novel stain
- Eliminate messy nail polish

745-1 Mycology Mounting Kit BLUE  
1 Fungitape + 12ml MycoPerm BLUE/pkg

745-9 Mycology Mounting Kit RED  
1 Fungitape + 12ml MycoPerm RED/pkg

752 MycoPerm RED - 12ml bottle/pkg


753 MycoPerm BLUE - 12ml bottle/pkg



### The Hines Method

*... an innovative technique*

The Mycology Laboratory at Hines VA Medical Center in suburban Chicago just loves their MycoPerm. Their innovative technique involves one drop of MycoPerm RED with acid fuchsin and one drop of MycoPerm BLUE with Lacto Phenol Cotton blue on each slide. As they scan the permanent mount they have a spectrum of color where fungal structures are sharper and even easier to visualize.



## Fungi-Tape

- Improve tape-preps for fungi slide mounts with specially formulated tape that does not disintegrate, fog or bubble
- Select colonies with easy to dispense, pliable adhesive tape - just form loop
- Professionalize process with strong tape, each roll contains approximately 360 coverslips

745 Fungi-Tape - 1 roll/pkg