

## Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g), Rev. 2012 and GHS Rev 03.

Printing date 05/15/2017

Reviewed on 05/15/2017

### \* 1 Identification

· **Product identifier**

· **Trade name: Helminth Ova and Larvae Formalin Suspensions**

· **Product number:**

| Cat# | Description                                    | SDL Prod ID |
|------|--|-------------|
| 510  | E. vermicularis Suspension – 1cc vial/pkg      | 00303       |
| 511  | Trichuris Suspension- 1cc vial/pkg             | 00304       |
| 512  | Strongyloides larvae Suspension - 1cc vial/pkg | 00305       |
| 515  | N. americanus Suspension - 1cc vial/pkg        | 00308       |
| 518  | H. nana Suspension - 1cc vial/pkg              | 00311       |
| 520  | S. japonicum Suspension - 1cc vial/pkg         | 00313       |
| 522  | Ascaris lumbricoides Suspension - 1cc vial/pkg | 00315       |

· **Details of the supplier of the safety data sheet**

· **Manufacturer/Supplier:**

Scientific Device Laboratory  
 411 Jarvis Avenue  
 Des Plaines, IL 60018  
 847-803-9495

· **Emergency telephone number:**

For Chemical Emergency Spill, Leak, Fire, Exposure, or Accident  
 Call CHEMTREC Day or Night  
 Within USA and Canada: 1-800-424-9300 or Outside USA and Canada: +1 703-527-3887 (collect calls accepted)

### \* 2 Hazard(s) identification

· **Classification of the substance or mixture**



GHS08 Health hazard

Carc. 2 H351 Suspected of causing cancer.



GHS07

Skin Irrit. 2 H315 Causes skin irritation.  
 Eye Irrit. 2A H319 Causes serious eye irritation.  
 Skin Sens. 1 H317 May cause an allergic skin reaction.

H227 Combustible liquid.

· **Label elements**

· **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).

· **Hazard pictograms**



GHS07 GHS08

**Safety Data Sheet (SDS)**

OSHA HazCom Standard 29 CFR 1910.1200(g), Rev. 2012 and GHS Rev 03.

Printing date 05/15/2017

Reviewed on 05/15/2017

**Trade name: Helminth Ova and Larvae formalin suspensions**

(Contd. of page 1)

- **Signal word** Warning
- **Hazard-determining components of labeling:**  
formaldehyde
- **Hazard statements**  
Combustible liquid.  
Causes skin irritation.  
Causes serious eye irritation.  
May cause an allergic skin reaction.  
Suspected of causing cancer.
- **Precautionary statements**  
Keep away from flames and hot surfaces. – No smoking.  
Avoid breathing dust/fume/gas/mist/vapors/spray  
Wear protective gloves / eye protection / face protection.  
Wear protective gloves.  
Wear eye protection / face protection.  
Wash thoroughly after handling.  
Contaminated work clothing must not be allowed out of the workplace.  
Obtain special instructions before use.  
Do not handle until all safety precautions have been read and understood.  
If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.  
Continue rinsing.  
Specific treatment (see on this label).  
Wash contaminated clothing before reuse.  
If exposed or concerned: Get medical advice/attention.  
If skin irritation or rash occurs: Get medical advice/attention.  
If eye irritation persists: Get medical advice/attention.  
In case of fire: Use for extinction: CO<sub>2</sub>, powder or water spray.  
If on skin: Wash with plenty of water.  
Take off contaminated clothing and wash it before reuse.  
Store locked up.  
Store in a well-ventilated place. Keep cool.  
Dispose of contents/container in accordance with local/regional/national/international regulations.

**\* 3 Composition/information on ingredients**

- **Chemical characterization: Mixtures**
- **Description:** Mixture of the substances listed below with nonhazardous additions.

(Contd. on page 3)

**Safety Data Sheet (SDS)**

OSHA HazCom Standard 29 CFR 1910.1200(g), Rev. 2012 and GHS Rev 03.

Printing date 05/15/2017

Reviewed on 05/15/2017

**Trade name: Helminth Ova and Larvae formalin suspensions**

(Contd. of page 1)

|                                |  |
|--------------------------------|--|
| <b>· Dangerous Components:</b> |  |
| 50-00-0                        | formaldehyde   |
|                                | ☠ Acute Tox. 3, H301; Acute Tox. 3, H311; Acute Tox. 3, H331; ☠ Carc. 2, H351; ☠ Skin Corr. 1B, H314; ⚠ Skin Sens. 1, H317 |

**4 First-aid measures****· Description of first aid measures**

- **After inhalation:** In case of unconsciousness, place patient securely on side position for transportation.
- **After skin contact:** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact:**  
Rinse opened eye for at least 15 minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:** If symptoms persist consult doctor.
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

**5 Fire-fighting measures**

- **Extinguishing media**
- **Suitable extinguishing agents:**  
CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Special hazards arising from the substance or mixture**  
Combustible liquid. Vapors can travel to a source of ignition and flash back.  
Explosive mixtures may occur at temperatures at or above flashpoint.
- **Advice for firefighters**
- **Protective equipment:**  
As in any fire, wear self-contained breathing apparatus pressure-demand (NIOSH approved or equivalent), and full protective gear to prevent contact with skin and eyes.

**\* 6 Accidental release measures**

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:**  
Dilute with plenty of water.  
Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**  
Absorb with liquid-binding material (ie. sand, diatomite, acid binders, universal binders, sawdust).  
Dispose contaminated material as waste according to section 13.  
Ensure adequate ventilation.  
Dispose of the collected material according to regulations.
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

**7 Handling and storage**

- **Precautions for safe handling**  
Ensure good ventilation/exhaustion at the workplace.  
Prevent formation of aerosols.
- **Information about protection against explosions and fires:** Protect from heat.

(Contd. on page 4)

## Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g), Rev. 2012 and GHS Rev 03.

Printing date 05/15/2017

Reviewed on 05/15/2017

**Trade name: Helminth Ova and Larvae formalin suspensions**

(Contd. of page 1)

- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:**  
 Keep receptacle tightly sealed.  
 Protect from heat and direct sunlight.
- **Specific end use(s)** No further relevant information available.

### 8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see section 7.
- **Control parameters**

· **Components with occupational exposure limits:**

**50-00-0 formaldehyde**

|     |   |
|-----|---|
| PEL | Short-term value: 2 ppm<br>Long-term value: 0.75 ppm<br>see 29 CFR 1910.1048(c)                 |
| REL | Long-term value: 0.016 ppm<br>Ceiling limit value: 0.1* ppm<br>*15-min; See Pocket Guide App. A |
| TLV | Ceiling limit value: 0.37 mg/m <sup>3</sup> , 0.3 ppm<br>(SEN) NIC-DSEN; RSEN                   |

- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**  
 Keep away from foodstuffs, beverages and feed.  
 Immediately remove all soiled and contaminated clothing.  
 Wash hands before breaks and at the end of work.  
 Avoid contact with the eyes and skin.
- **Breathing equipment:** Not required.
- **Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.  
 Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Select glove material based on penetration times, rates of diffusion and degradation.

· **Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact break-through time has to be determined and observed by the manufacturer of the protective gloves.

· **Eye protection:**

Safety glasses

(Contd. on page 5)

**Safety Data Sheet (SDS)**

OSHA HazCom Standard 29 CFR 1910.1200(g), Rev. 2012 and GHS Rev 03.

Printing date 05/15/2017

Reviewed on 05/15/2017

**Trade name: Helminth Ova and Larvae formalin suspensions**

(Contd. of page 1)

Tightly sealed goggles

**\* 9 Physical and chemical properties****· Information on basic physical and chemical properties****· General Information****· Appearance:**

**Form:** Liquid  
**Color:** Clear  
**Odor:** Pungent  
**Odor threshold:** Not determined.

**· pH-value:** Not determined.

**· Change in condition**

**Melting point/Melting range:** Not determined.  
**Boiling point/Boiling range:** -21 °C (-6 °F)

**· Flash point:** 85 °C (185 °F)

**· Flammability (solid, gaseous):** Not applicable.

**· Ignition temperature:** 300 °C (572 °F)

**· Decomposition temperature:** Not determined.

**· Auto igniting:** Product is not selfigniting.

**· Danger of explosion:** Not determined.

**· Explosion limits:**

**Lower:** 7.0 Vol %  
**Upper:** 73.0 Vol %

**· Vapor pressure @ 20 °C (68 °F):** 1.7 hPa (1 mm Hg)  
19 mmHg

**· Density @ 20 °C (68 °F):** 1.109 g/cm<sup>3</sup> (9.255 lbs/gal)

**· Relative density** Not determined.

**· Vapour density** 1.1 (AIR=1)

**· Evaporation rate** Not determined.

**· Solubility in / Miscibility with**

**Water:** Fully miscible.

**· Partition coefficient (n-octanol/water):** Not determined.

**· Viscosity:**

**Dynamic:** Not determined.  
**Kinematic:** Not determined.

**· Solvent content:**

**Organic solvents:** 10.0 %

**VOC content:** 10.0 %

**· Other information** No further relevant information available.

**10 Stability and reactivity**

**· Reactivity** No further relevant information available.

**· Chemical stability** Stable under normal conditions.

**· Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.

(Contd. on page 6)

## Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g), Rev. 2012 and GHS Rev 03.

Printing date 05/15/2017

Reviewed on 05/15/2017

**Trade name: Helminth Ova and Larvae formalin suspensions**

(Contd. of page 1)

- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

### \* 11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**

- **LD/LC50 values that are relevant for classification:**

**50-00-0 formaldehyde**

Oral LD50 >200 mg/kg (rat)

- **Primary irritant effect:**

- **on the skin:**

Irritant to skin and mucous membranes.  
 May cause an allergic skin reaction.

- **on the eye:**

Irritating effect.  
 Causes serious eye irritation.

- **Sensitization:** Sensitization possible through skin contact.

- **Additional toxicological information:**

The product shows the following dangers according to internally approved calculation methods for preparations:  
 Irritant

- **Carcinogenic categories**

- **IARC (International Agency for Research on Cancer)**

Group 1 - Carcinogenic to humans  
 Group 2A - Probably carcinogenic to humans  
 Group 2B - Possibly carcinogenic to humans  
 Group 3 - Not classifiable as to its carcinogenicity to humans  
 Group 4 - Probably not carcinogenic to humans

|                      |   |
|----------------------|---|
| 50-00-0 formaldehyde | 1 |
|----------------------|---|

- **NTP (National Toxicology Program)**

|                      |   |
|----------------------|---|
| 50-00-0 formaldehyde | K |
|----------------------|---|

- **OSHA-Ca (Occupational Safety & Health Administration)**

All ingredients are listed.

### 12 Ecological information

- **Toxicity**

- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.

- **Additional ecological information:**

- **General notes:**

Do not allow product to reach ground water, water course or sewage system.  
 Danger to drinking water if even small quantities leak into the ground.

- **Results of PBT and vPvB assessment**

- **PBT:** Not applicable.
- **vPvB:** Not applicable.

(Contd. on page 7)

## Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g), Rev. 2012 and GHS Rev 03.

Printing date 05/15/2017

Reviewed on 05/15/2017

**Trade name: Helminth Ova and Larvae formalin suspensions**

(Contd. of page 1)

· **Other adverse effects** No further relevant information available.

### \* 13 Disposal considerations

· **Waste treatment methods**

· **Recommendation:**

Must not be disposed of together with household garbage. Do not allow product to reach sewage system.

· **Uncleaned packagings:**

· **Recommendation:** Disposal must be made according to official regulations.

· **Recommended cleansing agent:** Water, if necessary with cleansing agents.

### \* 14 Transport information

· **UN-Number**

· **DOT, ADR, ADN, IMDG, IATA** Non-Regulated Material

· **UN proper shipping name**

· **DOT, ADR, ADN, IMDG, IATA** Non-Regulated Material

· **Transport hazard class(es)**

· **DOT, ADR, ADN, IMDG, IATA**

· **Class** Non-Regulated Material

· **Packing group**

· **DOT, ADR, IMDG, IATA** Non-Regulated Material

· **Environmental hazards:**

· **Marine pollutant:** No

· **Special precautions for user** Not applicable.

· **Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code** Not applicable.

· **UN "Model Regulation":** -

### \* 15 Regulatory information

· **Safety, health and environmental regulations/legislation specific for the substance or mixture**

· **Sara**

· **Section 355 (extremely hazardous substances):**

All ingredients are listed.

· **Section 313 (Specific toxic chemical listings):**

All ingredients are listed.

· **TSCA (Toxic Substances Control Act):**

All ingredients are listed.

· **Proposition 65**

· **Chemicals known to cause cancer:**

All ingredients are listed.

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients are listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients are listed.

· **Chemicals known to cause developmental toxicity:**

None of the ingredients are listed.

(Contd. on page 8)



## Safety Data Sheet (SDS)

OSHA HazCom Standard 29 CFR 1910.1200(g), Rev. 2012 and GHS Rev 03.

Printing date 05/15/2017

Reviewed on 05/15/2017

**Trade name: Helminth Ova and Larvae formalin suspensions**

(Contd. of page 1)

|   |              |    |
|---|--------------|----|
| · <b>EPA (Environmental Protection Agency)</b>                            |              |    |
| 50-00-0   | formaldehyde | B1 |
| · <b>TLV (Threshold Limit Value established by ACGIH)</b>                 |              |    |
| 50-00-0   | formaldehyde | A2 |
| · <b>NIOSH-Ca (National Institute for Occupational Safety and Health)</b> |              |    |
| None of the ingredients is listed.  |              |    |

· **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).

· **Hazard pictograms**



GHS07 GHS08

· **Signal word** Warning

· **Hazard-determining components of labeling:**

formaldehyde

· **Hazard statements**

Combustible liquid.

Causes skin irritation.

Causes serious eye irritation.

May cause an allergic skin reaction.

Suspected of causing cancer.

· **Precautionary statements**

Keep away from flames and hot surfaces. – No smoking.

Avoid breathing dust/fume/gas/mist/vapors/spray

Wear protective gloves / eye protection / face protection.

Wear protective gloves.

Wear eye protection / face protection.

Wash thoroughly after handling.

Contaminated work clothing must not be allowed out of the workplace.

Obtain special instructions before use.

Do not handle until all safety precautions have been read and understood.

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

Specific treatment (see on this label).

Wash contaminated clothing before reuse.

IF exposed or concerned: Get medical advice/attention.

If skin irritation or rash occurs: Get medical advice/attention.

If eye irritation persists: Get medical advice/attention.

In case of fire: Use for extinction: CO<sub>2</sub>, powder or water spray.

If on skin: Wash with plenty of water.

Take off contaminated clothing and wash it before reuse.

Store locked up.

Store in a well-ventilated place. Keep cool.

Dispose of contents/container in accordance with local/regional/national/international regulations.

(Contd. on page 9)



## Safety Data Sheet (SDS)





OSHA HazCom Standard 29 CFR 1910.1200(g), Rev. 2012 and GHS Rev 03.

Printing date 05/15/2017

Reviewed on 05/15/2017

**Trade name: Helminth Ova and Larvae formalin suspensions**

(Contd. of page 1)

|  |   |
|--|---|
| <b>· National regulations:</b><br>The product is subject to be labeled according with the prevailing version of the regulations on hazardous substances. |   |
| <b>· State Right to Know</b>   |   |
| 50-00-0  | formaldehyde<br> Acute Tox. 3, H301; Acute Tox. 3, H311; Acute Tox. 3, H331;  Carc. 2, H351;  Skin Corr. 1B, H314;  Skin Sens. 1, H317 |
| All ingredients are listed.  |   |
| <b>· Chemical safety assessment:</b> A Chemical Safety Assessment has not been carried out.  |   |

### 16 Other information

The information and recommendations in this safety data sheet are, to the best of our knowledge, accurate as of the date of issue. Nothing herein shall be deemed to create warranty, expressed or implied and shall not establish a legally valid contractual relationship. It is the responsibility of the user to determine applicability of this information and the suitability of the material or product for any particular purpose.

· **Date of preparation / last revision** 11/06/2014 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

ADN: The European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

VOC: Volatile Organic Compounds (USA, EU)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

: Flammable liquids, Hazard Category 4

Acute Tox. 3: Acute toxicity, Hazard Category 3

Skin Corr. 1B: Skin corrosion/irritation, Hazard Category 1B

Skin Irrit. 2: Skin corrosion/irritation, Hazard Category 2

Eye Irrit. 2A: Serious eye damage/eye irritation, Hazard Category 2A

Skin Sens. 1: Sensitisation - Skin, Hazard Category 1

Carc. 2: Carcinogenicity, Hazard Category 2

· **\* Data compared to the previous version altered.**

SDS created by MSDS Authoring Services [www.msdsauthoring.com](http://www.msdsauthoring.com) (877) 204-9106