



Transport temperature-sensitive material with safety and ease.

# Thermo-Tote Variable Temperature Incubators




## Transport Samples Reliably

- Reach desired temperature in only a few minutes
- Maintain temperatures from 25° to 50°C during transit
- Heat or cool from 2° to 60°C with ThermoTote Deluxe



## Operate in the Field or Lab

- Collect samples and immediately incubate on site
- Plug into standard outlet or car cigarette lighter
- Use rechargeable-battery model when no power source available
- Conduct special temperature studies in the lab



## Save Money


- Avoid lost samples and duplicated work due to temperature instability
- Cost-effective and affordable for labs of all types and sizes



### ACCESSORIES

- AC/DC Converter
- 220v/240v International Converter
- Cigarette Lighter Adapter
- Petri Dish Racks
- Multiple Female Plug

*For more information see the catalog on our website.*



*Transport water quality testing plates easily.*



*Transport inoculated vials and tubes safely.*



*Transport inoculated plates with the proper temperatures needed for maximum recovery.*

	<b>Thermo-Tote 12</b>	<b>Thermo-Tote 24</b>	<b>Thermo-Tote Deluxe</b>
TEMP RANGE	25°C to 50°C HEAT only	25°C to 50°C HEAT only	2°C to 60°C HEAT & COOL
POWER	12 volts	12 volts	12 volts
INTERNAL DIMENSIONS	7.75" x 10.75" x 7"	7.75" x 10.75" x 13.25"	15" x 10" x 12.5"
EXTERNAL DIMENSIONS	13" x 9.5" x 8.75"	13" x 9.5" x 15"	23.75" x 15.75" x 18.5"
WEIGHT	4 pounds	4.75 pounds	16 pounds
INCLUDES	Cigarette lighter power adapter plug	Cigarette lighter power adapter plug	Cigarette lighter power adapter plug

CATALOG #	DESCRIPTION	PACKAGING
7145	Thermo-Tote 12 Incubator adjustable w/ Display	Each
7146	Thermo-Tote 24 Incubator adjustable w/ Display	Each
7151	Thermo-Tote Deluxe Incubator	Each
7155 INTL	220v/240v Converter International	Each
7155 USA	AC/DC Converter 120V for standard U.S. plugs	Each
7162	Cigarette Lighter Adapter for small & medium ThermoTotes	Each
7160	Cigarette Lighter Adapter for large ThermoTote	Each
7152-01	Multiple Female Plug	Each
7159	Plastic Test Tube Racks	4 Racks/pkg
7157-12	Thermo-Tote 12 Petri Dish Rack (4 RACKS)	4 Racks/pkg
7157-24	Thermo-Tote 24 Petri Dish Rack (4 RACKS)	4 Racks/pkg
7156	Thermo-Tote Battery Guard	Each
7158	Thermo-Tote Shelves	2 Shelves/pkg